

MRI




Magnetic Resonance Imaging



DSA




Digital Subtraction Angiography



DSA Syringe

KL-8019-DSA-150-Mark


Range of Uses: Medrad Incorporated, USA
 ■ Mark V High Pressure Injector
 Standard Components: 1-150ml Syringe, 1-Pipette



DSA Syringe

KL-8019-DSA-150-Mark A


Range of Uses: Medrad Incorporated, USA
 ■ Mark VII High Pressure Injector
 Standard Components: 1-150ml Syringe, 1-Pipette

MRI Syringe

KL-8019-MRI-65/115-Mark

Range of Uses: Medrad Incorporated, USA
 ■ Spectris & Solaris MRI High Pressure Injector
 Standard Components: 1-65ml Syringe, 1-115ml Syringe, 1-2500mm Y (T) type Connecting Tube, 2-Pipettes



MRI Syringe

KL-8019-MRI-65/65-Mark

Range of Uses: Medrad Incorporated, USA
 ■ Spectris&Solaris MRI High Pressure Injector
 Standard Components: 2-65ml Syringes, 1-2500mm Y (T) -type Connecting Tube, 2-Pipettes



MRI Syringe

KL-8019-MRI-60/60-LF

Range of Uses: Liebel-Flarsheim Co.,LLC,USA
 ■ Optistar LE& Optistar Elite MRI High Pressure Injector
 Standard Components: 2-60ml Syringes, 1-2500mm Y (T) type Connecting Tube, 2-Pipettes



DSA Syringe

KL-8019-DSA-150-LF


Range of Uses: Liebel-Flarsheim Co., LLC, USA
 ■ Illumena High Pressure Injector
 Standard Components: 1-150ml Syringe, 1-Pipette



DSA Syringe

KL-8019-DSA-125-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ 120S DSA High Pressure Injector
 Standard Components: 1-125ml Syringe, 1-Pipette



DSA Syringe

KL-8019-DSA-200-MT


Range of Uses: Medtron Germany
 ■ Accutron HP-D DSA High Pressure Injector
 Standard Components: 1-200ml Syringe, 1-Pipette




MRI Syringe

KL-8019-MRI-50/50-NE


Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Sonic Shot 50 MRI High Pressure Injector
 Standard Components: 2-50ml Syringes, 1-2500mm Y (T) type Connecting Tube, 2-Pipettes



MRI Syringe

KL-8019-MRI-65/65-MT

Range of Uses: Medtron ,Germany
 ■ Accutron MRI High Pressure Injector
 Standard Components: 2-65ml Syringes, 1-2500mm Y (T) -type Connecting Tube, 2-Pipettes



MRI Syringe


KL-8019-MRI-100/100-EZ

Range of Uses: E-Z-EM Inc., USA
 ■ Bracco Acist MRI High Pressure Injector
 Standard Components: 2-100ml Syringes, 1-2500mm Y (T) -type Connecting Tube, 2-Pipettes




Accessory

Jiangsu Kangle Medical Devices Co.,Ltd.



Pipettes

J Quick Fill Tubing Type Spike Type



Tube System with Valves

24-hour technical support
 WhatsApp
 +61 430 464 796

一次性使用


高压造影注射器及附件

Disposable use High Pressure Angiographic Syringe and Accessories

Disposable use
**High Pressure Angiographic
Syringe and accessories**




**CT Dual Syringe
Computed Tomography**



CT (Dual Syringe)

KL-8019-CT-200/200-Mark B

Range of Uses: Medrad Incorporated, USA
 ■ Stellant -D CT High Pressure Injector
 Standard Components: 2-200ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-200/200-Mark C

Range of Uses: Imaxeon Pty Ltd, Australia
 ■ Imaxeon Salient High Pressure Injector
 Standard Components: 2-190ml syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



**CT Single Syringe
Computed Tomography**



CT (Single Syringe)

KL-8019-CT-200-Mark A

Range of Uses: Medrad Incorporated, USA
 ■ MCT&MCT Plus CT High Pressure Injector
 ■ Vistron CT High Pressure Injector
 ■ Envision High Pressure Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-Mark B

Range of Uses: Medrad Incorporated, USA
 ■ Stellant-S CT High Pressure Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-Mark C

Range of Uses: Imaxeon Pty Ltd, Australia
 ■ Imaxeon Salient High Pressure Injector
 Standard Components: 1-190ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-LF

Range of Uses: Liebel-Flarsheim Co., LLC, USA
 ■ CT9000& CT9000 ADV CT High Pressure
 Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Dual Syringe)

KL-8019-CT-200/200-LF

Range of Uses: Liebel-Flarsheim Co., LLC, USA
 ■ OptiVantage DH CT High Pressure Injector
 Standard Components: 2-200ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-200/200-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Dual Shot CT High Pressure Injector
 Standard Components: 2-200ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-100/100-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Dual Shot CT High Pressure Injector
 Standard Components: 2-100ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



CT (Single Syringe)

KL-8019-CT-100-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd,Japan
 ■ A-25& A-60 CT High Pressure Injector
 Standard Components: 1-100ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd,Japan
 ■ A-25& A-60 CT High Pressure Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-MT

Range of Uses: Medtronic Germany
 ■ Accutron CT High Pressure Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting tube, 1-Pipette



CT (Single Syringe)

KL-8019-CT-200-EZ


Range of Uses: E-Z-EM Inc., USA
 ■ Bracco Acist / Empower Pressure Injector
 Standard Components: 1-200ml Syringe,
 1-1500mm Connecting Tube, 1-Pipette



CT (Dual Syringe)

KL-8019-CT-100/200-NE


Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Dual Shot CT High Pressure Injector
 Standard Components: 1-100ml Syringe,
 1-200ml Syringe,1-1500mm Y (T) -type
 Connecting Tube, 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-50/100-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Dual Shot CT High Pressure Injector
 Standard Components: 1-50ml Syringe,
 1-100ml Syringe, 1-1500mm Y (T) -type
 Connecting Tube, 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-50/200-NE

Range of Uses: Nemoto Kyorindo Co.,Ltd, Japan
 ■ Dual Shot CT High Pressure Injector
 Standard Components: 1-50ml Syringe,
 1-200ml Syringe, 1-1500mm Y (T) -type
 Connecting Tube, 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-200/200-MT

Range of Uses: Medtronic Germany
 ■ Accutron CT -D High Pressure Injector
 Standard Components: 2-200ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes



CT (Dual Syringe)

KL-8019-CT-200/200-EZ

Range of Uses: E-Z-EM Inc. USA
 ■ Bracco Acist / Empower CT High
 Pressure Injector
 Standard Components: 2-200ml Syringes,
 1-1500mm Y (T) -type Connecting Tube,
 2-Pipettes